

Product datasheet for **HP205944**

CD51 (ITGAV) Human qPCR Primer Pair (NM_002210)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3685
Forward Sequence:	AGGAGAAGGTGCCTACGAAGCT
Reverse Sequence:	GCACAGGAAAGTCTTGCTAAGGC
ACCN:	NM_002210 , NM_002210.1 , NM_002210.3 , NM_002210.4 , BC047454 , BC072686 , BC126231 , BC136442 , BC144100 , BQ101180 , BQ129340
UniProt ID:	P06756
Synonyms:	CD51; MSK8; VNRA; VTNR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)