

Product datasheet for **HP205981**

KCNS1 Human qPCR Primer Pair (NM_002251)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3787
Forward Sequence:	TCTTCCGCGTACTCAAGTTGGC
Reverse Sequence:	GCCACACCAGAGAACTGACA
ACCN:	NM_002251 , NM_002251.1 , NM_002251.2 , NM_002251.3 , NM_002251.4 , BC075033 , BC075034 , BC132957 , BC132959 , BC144234 , NM_002251.5
UniProt ID:	Q96KK3
Synonyms:	Kv9.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)