

Product datasheet for **HP206000**

Cytokeratin 4 (KRT4) Human qPCR Primer Pair (NM_002272)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3851
Forward Sequence:	GCCGAGAATGACTTTGTGGTCC
Reverse Sequence:	CTCCGCATCATAGAGGACCTTC
ACCN:	<u>NM_002272</u> , <u>NM_002272.1</u> , <u>NM_002272.2</u> , <u>NM_002272.3</u> , <u>BC003630</u> , <u>BC012867</u> , <u>BC032132</u> , <u>BC042174</u>
UniProt ID:	<u>P19013</u>
Synonyms:	CK-4; CK4; CYK4; K4; WSN1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)