

Product datasheet for **HP206006**

KRTHA2 (KRT32) Human qPCR Primer Pair (NM_002278)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3882
Forward Sequence:	GGCAGAGAATGCCAGGATGGTT
Reverse Sequence:	GCACAGAGTGAGATCATCCAGG
ACCN:	<u>NM_002278</u> , <u>NM_002278.1</u> , <u>NM_002278.2</u> , <u>NM_002278.3</u> , <u>BC146286</u> , <u>BC148740</u>
UniProt ID:	<u>Q14532</u>
Synonyms:	HA2; hHa2; HKA2; KRTHA2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)