

Product datasheet for **HP206019**

Laminin gamma 1 (LAMC1) Human qPCR Primer Pair (NM_002293)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3915
Forward Sequence:	CTGTGAGGTCAACCACTTTGGG
Reverse Sequence:	AGCCTTCTCTGCATTACAGCG
ACCN:	<u>NM_002293</u> , <u>NM_002293.1</u> , <u>NM_002293.2</u> , <u>BC015586</u> , <u>BC035429</u> , <u>BC067813</u> , <u>BF507385</u> , <u>BI768612</u> , <u>NM_002293.4</u>
UniProt ID:	<u>P11047</u>
Synonyms:	LAMB2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)