

Product datasheet for **HP206022**

Plastin L (LCP1) Human qPCR Primer Pair (NM_002298)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3936
Forward Sequence:	GCTCTTGCTGAGGTGGGCTAAT
Reverse Sequence:	TCTCCTTTGGAGCCACCTGCT
ACCN:	<u>NM_002298</u> , <u>NM_002298.1</u> , <u>NM_002298.2</u> , <u>NM_002298.3</u> , <u>NM_002298.4</u> , <u>BC010271</u> , <u>BC010271.1</u> , <u>BC007673</u> , <u>BC015001</u> , <u>BC053839</u> , <u>BF893502</u>
UniProt ID:	<u>P13796</u>
Synonyms:	CP64; HEL-S-37; L-PLASTIN; LC64P; LPL; PLS2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)