

Product datasheet for **HP206048**

LTBR Human qPCR Primer Pair (NM_002342)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4055
Forward Sequence:	CAACGAGCACTGGAACCTACCTG
Reverse Sequence:	GGTCTTCGGTTTGCTTGTCAG
ACCN:	<u>NM_002342</u> , <u>NM_002342.1</u> , <u>NM_002342.2</u> , <u>BC026262</u> , <u>BC026262.1</u> , <u>BX648074</u> , <u>NM_002342.3</u>
UniProt ID:	<u>P36941</u>
Synonyms:	D12S370; LT-BETA-R; TNF-R-III; TNFCR; TNFR-RP; TNFR2-RP; TNFR3; TNFRSF3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)