

Product datasheet for **HP206055**

TACD2 (TACSTD2) Human qPCR Primer Pair (NM_002353)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4070
Forward Sequence:	GGACATCAAGGGCGAGTCTCTA
Reverse Sequence:	AGGCGCTTCATGGAGAACTTCG
ACCN:	<u>NM_002353</u> , <u>NM_002353.1</u> , <u>NM_002353.2</u> , <u>BC009409</u> , <u>BC009409.1</u> , <u>BM765888</u> , <u>BT007255</u> , <u>NM_002353.3</u>
UniProt ID:	<u>P09758</u>
Synonyms:	EGP-1; EGP1; GA733-1; GA7331; GP50; M1S1; TROP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)