

Product datasheet for **HP206058**

MARCKS Human qPCR Primer Pair (NM_002356)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 4082 |
| Forward Sequence: | CTCCTCGACTTCTTCGCCCAAG |
| Reverse Sequence: | TCTTGAAGGAGAAGCCGCTCAG |
| ACCN: | NM_002356 , NM_002356.1 , NM_002356.2 , NM_002356.3 , NM_002356.4 , NM_002356.5 , NM_002356.6 , BC089040 , BC089040.1 , BC013004 , BC028284 , BC039492 , BC068059 , BU678392 , NM_002356.7 |
| UniProt ID: | P29966 |
| Synonyms: | 80K-L; MACS; PKCSL; PRKCSL |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)