

Product datasheet for **HP206063**

MAG Human qPCR Primer Pair (NM_002361)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4099
Forward Sequence:	TGGTCGCCTTTGCCATCCTGAT
Reverse Sequence:	AGATGCGGAAGTCGCTGCTGAA
ACCN:	<u>NM_002361</u> , <u>NM_002361.1</u> , <u>NM_002361.2</u> , <u>NM_002361.3</u> , <u>BC053347</u> , <u>BC053347.1</u> , <u>BC034042</u> , <u>BC093045</u> , <u>BI755065</u> , <u>NM_002361.4</u>
UniProt ID:	<u>P20916</u>
Synonyms:	GMA; S-MAG; SIGLEC-4A; SIGLEC4A; SPG75
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)