

Product datasheet for **HP206073**

MARK3 Human qPCR Primer Pair (NM_002376)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4140
Forward Sequence:	CAGTCTCCTCACCACAAAGTGC
Reverse Sequence:	TGCTGGTCTGACTCCTTTTCGG
ACCN:	<u>NM_002376</u> , <u>NM_002376.1</u> , <u>NM_002376.2</u> , <u>NM_002376.5</u> , <u>NM_002376.6</u> , <u>BC024773</u> , <u>BC024773.1</u> , <u>BX161395</u> , <u>BX248250</u>
UniProt ID:	<u>P27448</u>
Synonyms:	CTAK1; KP78; Par-1a; PAR1A; VIPB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)