

## Product datasheet for **HP206079**

### MAX Human qPCR Primer Pair (NM\_002382)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4149
Forward Sequence:	TTGACGACCTCAAGCGGCAGAA
Reverse Sequence:	TTGGCGTTGGTGTAGAGGCTGT
ACCN:	<u><a href="#">BC003525</a></u> , <u><a href="#">BC004516</a></u> , <u><a href="#">NM_002382</a></u> , <u><a href="#">NM_002382.1</a></u> , <u><a href="#">NM_002382.2</a></u> , <u><a href="#">NM_002382.3</a></u> , <u><a href="#">NM_002382.4</a></u> , <u><a href="#">BC003525.1</a></u> , <u><a href="#">BC004516.2</a></u> , <u><a href="#">BC013669</a></u> , <u><a href="#">BC036092</a></u> , <u><a href="#">BC036092.1</a></u> , <u><a href="#">BC025685</a></u> , <u><a href="#">BC025685.1</a></u> , <u><a href="#">BC027924</a></u> , <u><a href="#">BF038005</a></u> , <u><a href="#">BF062526</a></u> , <u><a href="#">BG177865</a></u> , <u><a href="#">NM_002382.5</a></u>
UniProt ID:	<u><a href="#">P61244</a></u>
Synonyms:	bHLHd4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)