

Product datasheet for **HP206098**

MGAT5 Human qPCR Primer Pair (NM_002410)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4249
Forward Sequence:	AGGCAGAACCAGTCCCTTGTGT
Reverse Sequence:	CTTTGTGCTGGAGCCATAAACAG
ACCN:	NM_002410 , NM_002410.1 , NM_002410.2 , NM_002410.3 , NM_002410.4 , BC041917 , BC156313 , BC157048 , BU191961 , BX647087 , NM_002410.5
UniProt ID:	Q09328
Synonyms:	glcNAc-T V; GNT-V; GNT-VA; MGAT5A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)