

Product datasheet for **HP206102**

MIF Human qPCR Primer Pair (NM_002415)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4282
Forward Sequence:	AGAACCGCTCCTACAGCAAGCT
Reverse Sequence:	GGAGTTGTTCCAGCCACATTG
ACCN:	<u>NM_002415</u> , <u>NM_002415.1</u> , <u>BC022414</u> , <u>BC022414.1</u> , <u>BC000447</u> , <u>BC007376</u> , <u>BC007676</u> , <u>BC008914</u> , <u>BC013976</u> , <u>BC051440</u> , <u>BC053376</u> , <u>BT007148</u> , <u>NM_002415.2</u>
UniProt ID:	<u>P14174</u>
Synonyms:	GIF; GLIF; MMIF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)