

Product datasheet for **HP206120**

MPV17 Human qPCR Primer Pair (NM_002437)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4358
Forward Sequence:	GCAAGTTAGCCAACTTCTACCTGG
Reverse Sequence:	GCAAGGTGGAACGATGGAGTG
ACCN:	NM_002437 , NM_002437.1 , NM_002437.2 , NM_002437.3 , NM_002437.4 , BC001115 , BC001115.2 , BC016289 , BG714998 , BI602840 , NM_002437.5
UniProt ID:	P39210
Synonyms:	CMT2EE; MTDPS6; SYM1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)