

Product datasheet for **HP206121**

Mannose Receptor (MRC1) Human qPCR Primer Pair (NM_002438)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4360
Forward Sequence:	AGCCAACACCAGCTCCTCAAGA
Reverse Sequence:	CAAACGCTCGCGCATTGTCCA
ACCN:	<u>NM_002438</u> , <u>NM_002438.1</u> , <u>NM_002438.2</u> , <u>NM_002438.3</u> , <u>BC016003</u> , <u>BC142642</u> , <u>BC146838</u> , <u>BU198859</u>
UniProt ID:	<u>P22897</u>
Synonyms:	bA541119.1; CD206; CLEC13D; CLEC13DL; hMR; MMR; MRC1L1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)