

## Product datasheet for **HP206158**

### Myosin Phosphatase 2 (PPP1R12B) Human qPCR Primer Pair (NM\_002481)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4660
Forward Sequence:	CCATCTATCGCTCCTCTTCAAGC
Reverse Sequence:	CACTTGTGCTGTTGAGCTTCCG
ACCN:	<u><a href="#">NM_002481</a></u> , <u><a href="#">NM_002481.1</a></u> , <u><a href="#">NM_002481.2</a></u> , <u><a href="#">NM_002481.3</a></u> , <u><a href="#">BC110907</a></u> , <u><a href="#">BC034430</a></u> , <u><a href="#">BC071166</a></u> , <u><a href="#">BC144699</a></u> , <u><a href="#">BX501884</a></u> , <u><a href="#">NM_002481.4</a></u>
UniProt ID:	<u><a href="#">O60237</a></u>
Synonyms:	MYPT2; PP1bp55
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)