

## Product datasheet for **HP206164**

### NDUFA2 Human qPCR Primer Pair (NM\_002488)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4695
Forward Sequence:	TACGTGGAGCTGAAGAAGGCCGA
Reverse Sequence:	ACGTTCTCCAGGGCTCTGGTTA
ACCN:	<a href="#">NM_002488</a> , <a href="#">NM_002488.1</a> , <a href="#">NM_002488.2</a> , <a href="#">NM_002488.3</a> , <a href="#">NM_002488.4</a> , <a href="#">BC003674</a> , <a href="#">BC003674.1</a> , <a href="#">BM706616</a> , <a href="#">BQ054127</a> , <a href="#">NM_002488.5</a>
UniProt ID:	<a href="#">O43678</a>
Synonyms:	B8; CD14; CIB8; MC1DN13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)