

Product datasheet for **HP206169**

NDUFB6 Human qPCR Primer Pair (NM_002493)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4712
Forward Sequence:	TTCACTCATGTA CTGTACCTGTC
Reverse Sequence:	TCTGGACTTCTTTCAACTATGCC
ACCN:	NM_002493 , NM_002493.1 , NM_002493.2 , NM_002493.4 , BC009801 , BC016708 , BC029247 , BG220924 , BI669356 , BP208846 , BU952486 , NM_002493.5
UniProt ID:	O95139
Synonyms:	B17; CI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)