

Product datasheet for **HP206171**

NDUFS4 Human qPCR Primer Pair (NM_002495)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4724
Forward Sequence:	GAGATTGGCACAGGACCAGACT
Reverse Sequence:	CCAGACTGCATGTTATTGCGAGC
ACCN:	<u>NM_002495</u> , <u>NM_002495.1</u> , <u>NM_002495.2</u> , <u>NM_002495.3</u> , <u>BC005270</u> , <u>BC005270.1</u> , <u>BU145342</u> , <u>BU954271</u> , <u>NM_002495.4</u>
UniProt ID:	<u>O43181</u>
Synonyms:	AQDQ; CI-18; CI-18 kDa; CI-AQDQ; MC1DN1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)