

## Product datasheet for **HP206176**

### NEUROD1 Human qPCR Primer Pair (NM\_002500)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4760
Forward Sequence:	GGTGCCTTGCTATTCTAAGACGC
Reverse Sequence:	GCAAAGCGTCTGAACGAAGGAG
ACCN:	<u><a href="#">NM_002500</a></u> , <u><a href="#">NM_002500.1</a></u> , <u><a href="#">NM_002500.2</a></u> , <u><a href="#">NM_002500.3</a></u> , <u><a href="#">NM_002500.4</a></u> , <u><a href="#">BC009046</a></u> , <u><a href="#">BC009046.1</a></u> , <u><a href="#">BM505034</a></u> , <u><a href="#">BT019731</a></u>
UniProt ID:	<u><a href="#">Q13562</a></u>
Synonyms:	BETA2; BHF-1; bHLHa3; MODY6; NEUROD; T2D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)