

Product datasheet for **HP206183**

Nkx2.2 (NKX2-2) Human qPCR Primer Pair (NM_002509)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4821
Forward Sequence:	CCTTCTACGACAGCAGCGACAA
Reverse Sequence:	ACTTGAGCTTGAGTCCTGAGG
ACCN:	NM_002509 , NM_002509.1 , NM_002509.2 , NM_002509.3 , BC075093 , BC075093.2 , BC075092 , BM314196 , BQ632044 , BU070383 , NM_002509.4
UniProt ID:	O95096
Synonyms:	NKX2.2; NKX2B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)