

Product datasheet for **HP206190**

NPAS2 Human qPCR Primer Pair (NM_002518)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4862
Forward Sequence:	AAGTGCTGGGAACCTCAGGCTA
Reverse Sequence:	AGCCAGATCCACTGCTGACCTT
ACCN:	NM_002518 , NM_002518.1 , NM_002518.2 , NM_002518.3 , BC051351 , BC051351.1 , BC028107 , BC035412 , BC072383 , BM724652 , BX640699 , NM_002518.4
UniProt ID:	Q99743
Synonyms:	bHLHe9; MOP4; PASD4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)