

Product datasheet for **HP206195**

NRAS Human qPCR Primer Pair (NM_002524)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4893
Forward Sequence:	GAAACCTCAGCCAAGACCAGAC
Reverse Sequence:	GGCAATCCCATACAACCCCTGAG
ACCN:	<u>NM_002524</u> , <u>NM_002524.1</u> , <u>NM_002524.2</u> , <u>NM_002524.3</u> , <u>NM_002524.4</u> , <u>BC005219</u> , <u>BC005219.1</u> , <u>BC013214</u> , <u>BT019734</u> , <u>BT019735</u> , <u>NM_002524.5</u>
UniProt ID:	<u>P01111</u>
Synonyms:	ALPS4; CMNS; N-ras; NCMS; NRAS1; NS6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)