

## Product datasheet for **HP206206**

### LOX 1 (OLR1) Human qPCR Primer Pair (NM\_002543)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4973
Forward Sequence:	GAAACCCTTGCTCGGAAGCTGA
Reverse Sequence:	CAGATCCAGTCTTGCGGACAAG
ACCN:	<u><a href="#">NM_002543</a></u> , <u><a href="#">NM_002543.1</a></u> , <u><a href="#">NM_002543.2</a></u> , <u><a href="#">NM_002543.3</a></u> , <u><a href="#">BC022295</a></u> , <u><a href="#">BC055085</a></u> , <u><a href="#">BP219342</a></u> , <u><a href="#">BQ028409</a></u> , <u><a href="#">NM_002543.4</a></u>
UniProt ID:	<u><a href="#">P78380</a></u>
Synonyms:	CLEC8A; LOX1; LOXIN; SCARE1; SLOX1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)