

## Product datasheet for **HP206208**

### Osteoprotegerin (TNFRSF11B) Human qPCR Primer Pair (NM\_002546)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4982
Forward Sequence:	GGTCTCCTGCTAACTCAGAAAGG
Reverse Sequence:	CAGCAAACCTGAAGAATGCCTCC
ACCN:	<a href="#">BC030155</a> , <a href="#">NM_002546</a> , <a href="#">NM_002546.1</a> , <a href="#">NM_002546.2</a> , <a href="#">NM_002546.3</a> , <a href="#">BC011023</a> , <a href="#">BM748946</a> , <a href="#">BQ000897</a> , <a href="#">NM_002546.4</a>
UniProt ID:	<a href="#">O00300</a>
Synonyms:	OCIF; OPG; PDB5; TR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)