

Product datasheet for **HP206214**

ORC5L (ORC5) Human qPCR Primer Pair (NM_002553)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5001
Forward Sequence:	GTATCGTGGACAGCAGAGTTGC
Reverse Sequence:	CCATCAAGCTGATCGTCATGGC
ACCN:	NM_002553 , NM_002553.1 , NM_002553.2 , NM_002553.3 , BC023652 , BC023652.2 , BF678849 , NM_002553.4
UniProt ID:	O43913
Synonyms:	ORC5L; ORC5P; ORC5T; PPP1R117
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)