

Product datasheet for **HP206220**

P2RX4 Human qPCR Primer Pair (NM_002560)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5025
Forward Sequence:	GTGGCGGATTATGTGATACCAGC
Reverse Sequence:	CACACAGTGGTCGCATCTGGAA
ACCN:	<u>NM_002560</u> , <u>NM_002560.1</u> , <u>NM_002560.2</u> , <u>BC033826</u> , <u>BC033826.1</u> , <u>BC017458</u> , <u>BC035845</u> , <u>BT019738</u> , <u>BT019739</u> , <u>BX402887</u> , <u>NM_002560.3</u>
UniProt ID:	<u>Q99571</u>
Synonyms:	P2X4; P2X4R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)