

Product datasheet for **HP206240**

PARN Human qPCR Primer Pair (NM_002582)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5073
Forward Sequence:	CCTTGCGGAATTGGAAAAGCGG
Reverse Sequence:	CTCGTGGAGTTGTCAGAGGCT
ACCN:	NM_002582 , NM_002582.1 , NM_002582.2 , NM_002582.3 , BC050029 , BC050029.1 , BC044763 , BF435512 , NM_002582.4
UniProt ID:	O95453
Synonyms:	DAN; DKCB6; PFBMFT4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)