

Product datasheet for **HP206244**

PBX2 Human qPCR Primer Pair (NM_002586)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5089
Forward Sequence:	AAGCCTGCTCTCTTAGCGTCC
Reverse Sequence:	CTCTGCCAGAAGCATGTTGTCC
ACCN:	NM_002586 , NM_002586.1 , NM_002586.2 , NM_002586.3 , NM_002586.4 , BC082261 , BC082261.1 , BC003111 , BQ182713 , BT009823 , NM_002586.5
UniProt ID:	P40425
Synonyms:	G17; HOX12; PBX2MHC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)