

Product datasheet for **HP206246**

PCDHGC3 Human qPCR Primer Pair (NM_002588)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5098
Forward Sequence:	AACACGCTGCGGAGCTGTGATC
Reverse Sequence:	GCCTGAGAGAAACGCCAGTCC
ACCN:	NM_002588 , NM_002588.1 , NM_002588.2 , NM_002588.3 , BC019299 , BC004321 , BC006439 , BC026218 , BC027451 , BC027470 , BC063437 , BG765565 , NM_002588.4
UniProt ID:	Q9UN70
Synonyms:	PC43; PCDH-GAMMA-C3; PCDH2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)