

Product datasheet for **HP206252**

Proprotein Convertase 2 (PCSK2) Human qPCR Primer Pair (NM_002594)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5126
Forward Sequence:	CTTGGCTTCCACCTTCAGCAAC
Reverse Sequence:	GAATGCCTCAGAGTGCAGTTGC
ACCN:	<u>NM_002594</u> , <u>NM_002594.1</u> , <u>NM_002594.2</u> , <u>NM_002594.3</u> , <u>NM_002594.4</u> , <u>BC001905</u> , <u>BC005815</u> , <u>BC040546</u> , <u>BM893758</u> , <u>BP384064</u> , <u>NM_002594.5</u>
UniProt ID:	<u>P16519</u>
Synonyms:	NEC-2; NEC 2; NEC2; PC2; SPC2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)