

Product datasheet for **HP206255**

PHOS (PDC) Human qPCR Primer Pair (NM_002597)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5132
Forward Sequence:	CAGGCAAATGTCTTCTCCTCAGA
Reverse Sequence:	ACGGTATTTACGAAGGCAGTTTTTC
ACCN:	NM_002597 , NM_002597.1 , NM_002597.2 , NM_002597.3 , NM_002597.4 , BC093839 , BC112134 , BC148616 , BC153182 , BM674109 , BM715757 , BQ636049 , NM_002597.5
UniProt ID:	P20941
Synonyms:	MEKA; PHD; PhLOP; PhLP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)