

Product datasheet for **HP206270**

PEDF (SERPINF1) Human qPCR Primer Pair (NM_002615)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5176
Forward Sequence:	TGAAGCGAAGTCACCAAGTCC
Reverse Sequence:	CCATCCTCGTTCCACTCAAAGC
ACCN:	<u>NM_002615</u> , <u>NM_002615.1</u> , <u>NM_002615.2</u> , <u>NM_002615.3</u> , <u>NM_002615.4</u> , <u>NM_002615.5</u> , <u>NM_002615.6</u> , <u>BC013984</u> , <u>BC013984.2</u> , <u>BC000522</u> , <u>BM981180</u> , <u>BQ773314</u> , <u>BT007222</u> , <u>BU154385</u>
UniProt ID:	<u>P36955</u>
Synonyms:	EPC-1; OI6; OI12; PEDF; PIG35
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)