

Product datasheet for **HP206274**

PF4 Human qPCR Primer Pair (NM_002619)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5196
Forward Sequence:	TCCTGCCACTTGTGGTCGCCT
Reverse Sequence:	CCTTGATCACCTCCAGGCTGG
ACCN:	<u>NM_002619</u> , <u>NM_002619.1</u> , <u>NM_002619.2</u> , <u>NM_002619.3</u> , <u>BC112093</u> , <u>BC093965</u> , <u>NM_002619.4</u>
UniProt ID:	<u>P02776</u>
Synonyms:	CXCL4; PF-4; SCYB4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)