

Product datasheet for **HP206281**

PFKP Human qPCR Primer Pair (NM_002627)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5214
Forward Sequence:	AGGCAGTCATCGCCTTGCTAGA
Reverse Sequence:	ATCGCCTTCTGCACATCCTGAG
ACCN:	NM_002627 , NM_002627.1 , NM_002627.2 , NM_002627.3 , NM_002627.4 , BC002536 , BC002536.1 , BC021698 , BC029138 , BC043156 , BU737488
UniProt ID:	Q01813
Synonyms:	ATP-PFK; PFK-C; PFK-P; PFKF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)