

Product datasheet for **HP206284**

Pepsinogen II (PGC) Human qPCR Primer Pair (NM_002630)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5225
Forward Sequence:	ACCTACTCCACCAATGGGCAGA
Reverse Sequence:	TCACTCAAGCCGAACTCCTGGT
ACCN:	<u>NM_002630</u> , <u>NM_002630.1</u> , <u>NM_002630.2</u> , <u>NM_002630.3</u> , <u>BC073740</u> , <u>BC073740.1</u> , <u>BC042578</u> , <u>NM_002630.4</u>
UniProt ID:	<u>P20142</u>
Synonyms:	PEPC; PGII
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)