

Product datasheet for **HP206287**

PGM1 Human qPCR Primer Pair (NM_002633)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5236
Forward Sequence:	TGATGGACGCGAGCAAAGTGTGTC
Reverse Sequence:	ATGTCCTCCACTCTGCTTGC
ACCN:	<u>NM_002633</u> , <u>NM_002633.1</u> , <u>NM_002633.2</u> , <u>BC019920</u> , <u>BC019920.1</u> , <u>BC001756</u> , <u>BC067763</u> , <u>BC090856</u> , <u>BT006961</u> , <u>NM_002633.3</u>
UniProt ID:	<u>P36871</u>
Synonyms:	CDG1T; GSD14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)