

Product datasheet for **HP206304**

GCDFP 15 (PIP) Human qPCR Primer Pair (NM_002652)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5304
Forward Sequence:	CCTGCCTATGTGACGACAATCC
Reverse Sequence:	TCAGGGCAGATGCCTAATTCCC
ACCN:	BC010950 , NM_002652 , NM_002652.1 , NM_002652.2 , BC010950.1 , BC010951 , NM_002652.3
UniProt ID:	P12273
Synonyms:	GCDFP-15; GCDFP15; GPIP4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)