

Product datasheet for **HP206308**

uPA (PLAU) Human qPCR Primer Pair (NM_002658)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5328
Forward Sequence:	GGCTTAACCTCCAACACGCAAGG
Reverse Sequence:	CCTCCTTGGAACGGATCTTCAG
ACCN:	<u>NM_002658</u> , <u>NM_002658.1</u> , <u>NM_002658.2</u> , <u>NM_002658.3</u> , <u>NM_002658.4</u> , <u>BC013575</u> , <u>BC013575.2</u> , <u>BI759308</u> , <u>BI833559</u> , <u>BQ920793</u> , <u>BT007391</u> , <u>BX360462</u> , <u>NM_002658.5</u>
UniProt ID:	<u>P00749</u>
Synonyms:	ATF; BDPLT5; QPD; u-PA; UPA; URK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)