

Product datasheet for **HP206309**

Phospholipase C gamma 1 (PLCG1) Human qPCR Primer Pair (NM_002660)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5335
Forward Sequence:	CATCACGCACTACCAGCAGGTG
Reverse Sequence:	GACGCGCATTAGCATGTGCTCA
ACCN:	<u>NM_002660</u> , <u>NM_002660.1</u> , <u>NM_002660.2</u> , <u>BC136466</u> , <u>BC144136</u> , <u>BX415739</u> , <u>BX434193</u> , <u>NM_002660.3</u>
UniProt ID:	<u>P19174</u>
Synonyms:	NCKAP3; PLC-II; PLC1; PLC148; PLCgamma1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)