

Product datasheet for **HP206310**

PLCG 2 (PLCG2) Human qPCR Primer Pair (NM_002661)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5336
Forward Sequence:	AGCGACTCCTATGCCATCACCT
Reverse Sequence:	AACTGACGAGCTCCACCAGACT
ACCN:	<u>NM_002661</u> , <u>NM_002661.1</u> , <u>NM_002661.2</u> , <u>NM_002661.3</u> , <u>NM_002661.4</u> , <u>BC014561</u> , <u>BC014561.2</u> , <u>BC007565</u> , <u>BC011772</u> , <u>BC018646</u> , <u>BC018983</u> , <u>BP351485</u> , <u>BX504887</u> , <u>NM_002661.5</u>
UniProt ID:	<u>P16885</u>
Synonyms:	APLAID; FCAS3; PLC-gamma-2; PLC-IV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)