

Product datasheet for **HP206318**

PLRG1 Human qPCR Primer Pair (NM_002669)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5356
Forward Sequence:	CGATGTGTTGGTAACCTGTAGTC
Reverse Sequence:	CAGCCTGACATCTCACTGTAGC
ACCN:	<u>NM_002669</u> , <u>NM_002669.1</u> , <u>NM_002669.2</u> , <u>NM_002669.3</u> , <u>BC106004</u> , <u>BC106004.1</u> , <u>BC020786</u> , <u>BC036046</u> , <u>BP314297</u> , <u>NM_002669.4</u>
UniProt ID:	<u>O43660</u>
Synonyms:	Cwc1; PRL1; PRP46; PRPF46; TANGO4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)