

Product datasheet for **HP206320**

Plexin B1 (PLXNB1) Human qPCR Primer Pair (NM_002673)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5364
Forward Sequence:	TCTCACCTGAATGGCTCCAAG
Reverse Sequence:	CTGGTCTCACACCGCAGTTGTT
ACCN:	NM_002673 , NM_002673.1 , NM_002673.2 , NM_002673.3 , NM_002673.4 , NM_002673.5 , BC066773 , BC128094 , BC129947 , BC146793
UniProt ID:	O43157
Synonyms:	PLEXIN-B1; PLXN5; SEP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)