

## Product datasheet for **HP206348**

### PPP1R3A Human qPCR Primer Pair (NM\_002711)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5506
Forward Sequence:	CCTAACAGACAAATAAAAGGCTGC
Reverse Sequence:	GAGAACAATGATTGTTGGGATATAG
ACCN:	<a href="#">NM_002711</a> , <a href="#">NM_002711.1</a> , <a href="#">NM_002711.2</a> , <a href="#">NM_002711.3</a> , <a href="#">BC126451</a> , <a href="#">BC017405</a> , <a href="#">BC056666</a> , <a href="#">BC070313</a> , <a href="#">BC111452</a> , <a href="#">BC126453</a> , <a href="#">BF576879</a>
UniProt ID:	<a href="#">Q16821</a>
Synonyms:	GM; PP1G; PPP1R3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)