

Product datasheet for **HP206364**

Major Basic Protein (PRG2) Human qPCR Primer Pair (NM_002728)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5553
Forward Sequence:	CCTGGTGAGAAGTCTTCAGACG
Reverse Sequence:	CGCTGACAGAACACTGGATTCG
ACCN:	NM_002728 , NM_002728.1 , NM_002728.2 , NM_002728.3 , NM_002728.4 , NM_002728.5 , BC005929 , BC005929.1 , BQ445845 , BX395670 , BX397850 , BX454677 , NM_002728.6
UniProt ID:	P13727
Synonyms:	BMPG; MBP; MBP1; proMBP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)