

Product datasheet for **HP206408**

PSG4 Human qPCR Primer Pair (NM_002780)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5672
Forward Sequence:	AGTGACCCAGTCACCCTGAATG
Reverse Sequence:	ATTGTGCCCGTGGGTTAGACTC
ACCN:	NM_002780 , NM_002780.1 , NM_002780.2 , NM_002780.3 , NM_002780.4 , BC063127 , BC063127.1 , BC008405 , BC013431 , BC035419 , BG006351 , BX403624 , NM_002780.5
UniProt ID:	Q00888
Synonyms:	PSBG-4; PSBG-9; PSG9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)