

Product datasheet for **HP206425**

PSMB9 Human qPCR Primer Pair (NM_002800)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 5698 |
| Forward Sequence: | CGAGAGGACTTGTCTGCACATC |
| Reverse Sequence: | CACCAATGGCAAAGGCTGTCG |
| ACCN: | NM_002800 , NM_002800.1 , NM_002800.2 , NM_002800.3 , NM_002800.4 , BC065513 , BC065513.1 , BF513693 , BI832878 , BM998687 , BX641100 , NM_002800.5 |
| UniProt ID: | P28065 |
| Synonyms: | beta1i; LMP2; PRAAS3; PSMB6i; RING12 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)