

Product datasheet for **HP206447**

PTK9 (TWF1) Human qPCR Primer Pair (NM_002822)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5756
Forward Sequence:	CCCAAGGATTCAGCTCGTTACC
Reverse Sequence:	GCTAGAATACAGCATCCGCTCTC
ACCN:	<u>NM_002822</u> , <u>NM_002822.1</u> , <u>NM_002822.2</u> , <u>NM_002822.3</u> , <u>NM_002822.4</u> , <u>BC043148</u> , <u>BC022344</u> , <u>BC022491</u> , <u>BC068548</u> , <u>BT019691</u>
UniProt ID:	<u>Q12792</u>
Synonyms:	A6; PTK9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)