

Product datasheet for **HP206452**

PTPH1 (PTPN3) Human qPCR Primer Pair (NM_002829)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5774
Forward Sequence:	GGACATCTCAGAACACACGCATG
Reverse Sequence:	GAAGTCAGCAAATGAGCGGACAG
ACCN:	NM_002829 , NM_002829.1 , NM_002829.2 , NM_002829.3 , BC143848 , BC126117 , BC033716 , BC063287 , BU624870 , BX648253 , BX648735 , NM_002829.4
UniProt ID:	P26045
Synonyms:	PTP-H1; PTPH1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)